


**IC-MPPE / Integrated
Computational Materials
Process and Product
Engineering**

Program: COMET – Competence
Centers for Excellent Technologies

Program line: COMET Center (K2)

Type of project: Strategic COMET
project



Experimental setup of a scanning electron microscope in electron backscatter diffraction mode during EBSD measurement at 1000 °C

NOVEL CHARACTERIZATION TECHNIQUES FOR MODEL VALIDATION AND PARAMETER DETERMINATION

HIGH TEMPERATURE ELECTRON BACKSCATTER DIFFRACTION AND COMPUTER-ASSISTED CREEP FRACTURE MECHANICAL TESTING

Together with partners, the Materials Center Leoben has developed a variety of material models in recent years to predict process-microstructure-property relationships for developing new materials, manufacturing processes or components. As models need to become more complex and accurate, characterization techniques must also evolve to validate these new models and determine reliable model parameters. In the following, two novel characterization techniques for model validation and parameter determination are presented as examples.

High temperature electron backscatter diffraction for validation of austenite grain-growth models

A new austenite grain-growth model was developed in order to predict the austenite grain size evolution

during the austenitization treatment of a high-alloyed steel. In order to validate the austenite grain-growth model and to determine required input parameters, high temperature electron backscatter diffraction (HT-EBSD) measurements were conducted. The *in situ* HT-EBSD observation during heating to austenitization temperature, illustrated in Fig.1, reveals novel phenomena for this steel alloy such as the austenite memory effect (i.e. the formation of new austenite grains during austenitization with the same crystal orientation, size and shape as the prior austenite grain – Fig.1b) and subsequent spontaneous recrystallization (i.e. recrystallization of austenite without previous deformation – Fig.1c).

SUCCESS STORY




Fig.1: EBSD imaging of a high-alloyed steel: a) Crystal orientation (inverse pole figure color coding) of martensitic microstructure at room temperature, (b) corresponding austenitic microstructure after heating to intercritical annealing temperature indicating the austenite memory effect, and (c) spontaneous recrystallization after further heating to the austenitization temperature above 1000 °C.

The HT-EBSD results enabled the determination of the initial grain size (after spontaneous recrystallization) for the austenite grain-growth model and improved the prediction of austenite grain size evolution during austenitization. The implementation of a spontaneous recrystallization model is ongoing research.

Computer-aided characterization of creep fracture mechanical properties at elevated temperatures

The determination of fracture mechanical properties at elevated temperatures for model validation is challenging due to large-scale plasticity and creep. Creep causes time-dependent specimen conditions (i.e. strain and stress) and complicates classical

elasto-plastic fracture mechanics approaches. Therefore, a new computer-aided testing technique was developed that includes the calculation of time-dependent specimen conditions (e.g. equivalent creep strain) using finite element modeling (FEM). The experimental setup is shown in Fig.2a at room temperature and in Fig. 2b at the testing temperature of 700°C under constant tension of 25 kN. The area for the corresponding FEM calculations is highlighted by a yellow rectangle. The results of the FEM calculation for the time-dependent equivalent creep strain are shown in Fig. 2c for 0.54 h and Fig. 2d for 20 h. The combination of experiment and simulation allows the characterization of fracture mechanical properties at temperatures up to 1000°C by correcting the experimental data according to the creep strain.




Fig. 2: Experimental setup of computer-aided fracture mechanical testing of a CT specimen according to ATSM E399 with a thickness of 25 mm made of an austenitic heat resistant steel: (a) testing at room temperature, (b) testing at 700°C with the highlighted region for the FEM calculation. The calculated equivalent creep strain is illustrated for 0.54 h (c) and 20 h (d) with blue regions for values larger than 10^{-5} .

Project coordination (Story)

Dr. Dominik Brandl
Junior Researcher Steel Engineering
Department Materials
T +43 (0) 3842 45922 - 62
dominik.brandl@mcl.at

Materials Center Leoben Forschung GmbH

Coordinator: COMET K2 Center IC-MPPE
Roseggerstrasse 12
8700 Leoben
T +43 (0) 3842 45922-0
mclburo@mcl.at
www.mcl.at

Project partner

- Erich Schmid Institute of Materials Science, Austrian Academy of Sciences, Austria
- University of Maribor, Slovenia

This success story was provided by the consortium leader/center management and by the mentioned project partners for the purpose of being published on the FFG website. Further information on COMET: www.ffg.at/comet